### **ARTEL**





# MVS

Multichannel Verification System


Increase liquid handling quality with easy and reliable performance verifications.



## Streamline volume verifications with Artel MVS® and ArtelWare™ Software

Without knowing the exact volume transferred at each step of an assay, the concentration of species in solution is unknown, potentially leading to unreliable results, delays, and cost overruns. Utilization of the Artel MVS for performance verification can prevent problems, bolster laboratory productivity, and provide confidence in the integrity of data.

# Three simple steps to volume verification with automatic reporting.



1. Dispense the QualAssure solution and Diluent into a Verification Plate and mix



2. Measure the absorbance values with the Plate Reader



3. Software displays results that indicate the performance of the liquid handler





#### Ideal for:

- Verifying the accuracy and precision of automated liquid handlers and multichannel pipettes
- Comparing instrument performance regardless of make, model or location
- Obtaining highly accurate results traceable to national and international standards
- Optimizing custom liquid classes on an automated liquid handler

### Verify volume transfer precision and accuracy and simplify regulatory compliance

#### **System Capabilities:**

- Measures performance over a wide volume range of 0.0001 to 350 µL for Automated Liquid Handlers and Multichannel Pipettes
- Allows for aqueous and non-aqueous volume transfers with specialized solutions such as DMSO, PCRMix, SerumSub, PlasProxy and more
- · Obtains volume statistics tip-by-tip and well-by-well
- Eliminates the need for rigorous environmental conditions
- Simplifies verification and calibration and is easy-to-use for technicians of any skill level

 Generates measurement results traceable to the International System of Units (SI) via standards maintained by NIST

#### **Key Applications:**

- Optimizing liquid classes
- Facilitating of scale-up and method transfer
- Developing, troubleshooting, and validating assays, methods, and SOPs
- · Ensuring reproducibility of assay results
- Site acceptance testing and establishing baseline performance of new equipment
- · Verifying performance before and after maintenance
- Assessing performance of other essential liquid handling equipment such as plate shakers, plate washers, and bulk dispensers

#### Compliance:

The MVS conforms to the Dual-Dye Ratiometric Photometric Method described in ISO 23783-2:2022, and the Ratiometric Photometric Method described in ISO/IWA 15:2015 for automated liquid handling systems. It can be used to fulfill test and calibration requirements of handheld multichannel pipettes according to ISO 8655-7:2022. The MVS can be used to fulfill test and calibration requirements of handheld or automated liquid handling systems according to CLSI QMS23:2019, ISO 17025:2017, cGMP, and cGLP. ArtelWare and Data Manager Software can be used in laboratories compliant to 21 CFR Part 11.

"The MVS is easy to use, facilitating a frequent performance monitoring schedule. It has made a big difference in our laboratories, both from an efficiency point of view, as well as from a confidence standpoint."





#### **Parts and Supplies**

Artel 800TSNB Plate Reader			
Size	(DxWxH): 41.9 x 38.1 x 17.8 cm		
Weight	9.97 kg		
Disply	Touch screen		
Light Source	Tungsten gas filled bulb		
Wavelength Selection	Metal oxide interference filters Center wavelengths: 520.2 nm (bandwidth at half-max = 6.2 nm) 730.5 nm (bandwidth at half-max = 10 nm)		
Additional Filters	405 nm, 450 nm, 490 nm		
COM port	USB		
Electrical Power Requirements	Voltage 90 – 260 VAC		
Frequency	50 – 60Hz		
Current	2 A maximum		
Warranty	1 year		

QualAssure Solution Types:		
Aqueous		
Volume Range	0.0001 to 350.0 μL	
Shelf life	24 months*	
Storage/Operation Temp	15-30 °C	
DMSO		
Volume Range	0.0001 to 9.999 μL	
Shelf Life	24 months*	
Storage/Operation Temp	19-30 °C	
PCRMix		
Volume Range	2.000 to 49.99 μL	
Shelf Life	24 months*	
Storage/Operation Temp	2-25 °C	
SerumSub and PlasProxy		
Volume Range	10.00 to 200.0 μL	
Shelf Life	12 months*	
Storage/Operation Temp	2-8 °C	

<sup>\*</sup>From date of manufacture

#### **Specifications**

96 well MVS Verification Plate\*\*\*

Time Requirements		< 5 minutes	
Tip Configurations		1, 2, 4, 6, 8, 12, and 96	
	Solution Range (µL)	(Uncertainty) Inaccuracy (%)**	(Imprecision) CV (%)**
Aqueous HV	350.0 to 200.1	1.40 to 1.44	0.15 to 0.16
Aqueous A	200.0 to 50.00	1.36 to 1.58	0.16 to 0.19
Aqueous B	49.99 to 10.00	1.54 to 1.90	0.15 to 0.20
Aqueous C	9.999 to 2.000	1.46 to 1.63	0.15 to 0.20
Aqueous D	1.999 to 1.000	1.38 to 1.61	0.16 to 0.19
Aqueous E	0.9999 to 0.1000	1.39 to 2.87	0.16 to 0.22
	0.0999 to 0.0001 (Non- traceable volume range)*	N/A	N/A
Operating Volume Range		15 - 30 °C	

#### 384 well MVS Verification Plate\*\*\*

304 Well WV3 Verification Flate			
Time Requirements	< 10 minutes		
Tip Configurations	1, 8, 12, 16, 24, 96, and 384		
	Solution Range (μL)	(Uncertainty) Inaccuracy (%)**	(Imprecision) CV (%)**
Aqueous HV	N/A	N/A	N/A
Aqueous A	55.00 to 10.00	1.87 to 2.08	0.36 to 0.39
Aqueous B	9.999 to 2.500	2.02 to 2.33	0.36 to 0.40
Aqueous C	2.499 to 0.500	2.04 to 2.35	0.36 to 0.40
Aqueous D	0.4999 to 0.3000	1.88 to 2.01	0.36 to 0.38
Aqueous E	0.2999 to 0.0300	1.89 to 3.23	0.36 to 0.44
	0.0299 to 0.0100	3.24 to 7.47	0.44 to 0.47
	0.0099 to 0.0001 (Non-traceable volume range)	N/A	N/A
Operating Volume Range		15 - 30 °C	

- \* Measurement of volumes outside of the traceable volume ranges are not traceable to the national and international standards and no declarations of relative inaccuracy and imprecision are made.
- \*\*\* Stated specifications apply when Artel MVS Plate Readers (ELx800 and 800TSNB) and MVS Verification Plates are used with Aqueous QualAssure and DMSO QualAssure.
- \*\*\* Requires 384-Well Plate Support add-on activation code.





#### Optimize your automation with QualAssure solutions that mimic your liquid types

- Saves development time
- Ensures correct volume delivery
- Decreases assay variability



25 Bradley Drive, Westbrook, Maine 04092 USA | 888-406-3463 | aicompanies.com Global Headquarters: Two Technology Way | Norwood, Massachusetts 02062, USA © 2023 Advanced Instruments | Artel, MVS, PCS, and QualAssure are trademarks of Advanced Instruments.